



How To Use Aerolizer™

CUE WORDS

1. **REMOVE** the cap.
2. **OPEN** by turning the mouthpiece.
3. **PLACE** capsule in base of Aerolizer™.
4. **CLOSE** by turning back the mouthpiece.
5. **HOLD** inhaler upright; **FIRMLY SQUEEZE** two blue buttons once; release.
6. **EXHALE**, then **SEAL LIPS** around the mouthpiece.
7. **INHALE** quickly and deeply.
8. **HOLD BREATH** for 5 – 10 seconds.
9. **EXHALE**.
10. **OPEN** inhaler to ensure capsule is empty.
11. **INHALE AGAIN** if not.