



How To Use A Diskhaler®

CUE WORDS

1. **REMOVE** the cover.
2. **ADVANCE** the disk to the next dose.
3. **HOLD** the Diskhaler, flat.
4. **LIFT** the cover panel to **PIERCE** the shell.
5. **CLOSE** the cover.
6. **EXHALE**, then **SEAL LIPS** around mouthpiece.
7. **INHALE QUICKLY AND DEEPLY**.
8. **HOLD BREATH** for 5 – 10 seconds.
9. **EXHALE**.